

Non-Shouldered Abutment System

Implant System	3.5mm Diameter		4.0mm Diameter			5.0mm Diameter			5.0mm Diameter			6.5mm Diameter		7.5mm Diameter																								
	0°	15°	0°	15°	25°	0°	15°	25°	0°	15°	25°	0°	15°	0°	15°																							
SHORT IMPLANTS 3.0mm Post																																						
	4.0 x 6.5 0° 260-340-001		4.0 x 6.5 15° 260-340-015		4.0 x 10 0° 260-340-101		4.0 x 10 15° 260-340-115		5.0 x 5.0 0° 260-350-050		5.0 x 5.0 15° 260-350-055		5.0 x 6.5 0° 260-350-001		5.0 x 6.5 15° 260-350-015		5.0 x 6.5 25° 260-350-025		5.0 x 10 0° 260-350-101		5.0 x 10 15° 260-350-115		5.0 x 12 0° 260-350-201		5.0 x 12 15° 260-350-215		6.5 x 5.0 0° 260-365-050		6.5 x 5.0 15° 260-365-055		6.5 x 6.5 0° 260-365-001		6.5 x 6.5 15° 260-365-015		7.5 x 8.0 0° 260-375-801		7.5 x 8.0 15° 260-375-815	
MAX 2.5™ MAXILLARY ANTERIOR IMPLANTS 2.5mm Post																																						
	4.0 x 6.5 0° 260-240-001		4.0 x 6.5 15° 260-240-015		4.0 x 6.5 25° 260-240-025		4.0 x 10 0° 260-240-101		4.0 x 10 15° 260-240-115		5.0 x 6.5 0° 260-250-001		5.0 x 6.5 15° 260-250-015		5.0 x 6.5 25° 260-250-025		5.0 x 10 0° 260-250-111		5.0 x 10 15° 260-250-115		6.5 x 5.0 0° 260-265-050		6.5 x 5.0 15° 260-265-055															
NARROW IMPLANTS 2.0mm Post																																						
	3.5 x 6.5 0° 260-135-001		3.5 x 6.5 15° 260-135-015		4.0 x 6.5 0° 260-140-002		4.0 x 6.5 15° 260-140-015		4.0 x 6.5 25° 260-140-025		4.0 x 10 0° 260-140-101		4.0 x 10 15° 260-140-115		5.0 x 5.0 0° 260-150-050		5.0 x 5.0 15° 260-150-055		5.0 x 6.5 0° 260-150-001		5.0 x 6.5 15° 260-150-015		5.0 x 6.5 25° 260-150-025		5.0 x 12 0° 260-150-201		5.0 x 12 15° 260-150-215											
RESTORATIVE/LABORATORY KIT																																						
	3.5 x 6.5mm 260-135-465		4.0 x 6.5mm 260-140-465		4.0 x 10mm 260-140-410		5.0 x 5.0mm 260-150-450		5.0 x 6.5mm 260-150-465		5.0 x 10mm 260-150-410		5.0 x 12mm 260-150-412		6.5 x 5.0mm 260-165-450		6.5 x 6.5mm 260-165-465		7.5 x 8.0mm 260-175-480																			
	3.5mm 260-135-165		4.0mm 260-140-165				5.0mm 260-150-165				5.0mm 260-150-165				6.5mm 260-165-165				7.5mm 260-175-165																			

- **Notes:**
- Snap-on sleeves are only specific for abutment diameter.
 - Abutment height is not a criterion for proper selection of snap-on sleeves.
 - Transfer dies correspond to exact diameter and height of abutment placed.
 - Because of machining tolerances, sleeves may not reach the height of contour for some angled abutments.